



Large Letter Mailmark barcode quality and item design check

1. Provide 10 fully addressed, filled and sealed mail pieces with the appropriate indicia and Mailmark barcode. Items will be retained for future reference
2. Provide details of the type of wrapping used for the samples sent.
Please provide a separate sheet per type to enable us to differentiate during testing.

Wrap type	Quantity	Haze	Micron thickness
Polyethylene [PE]			
Polypropylene [PP]			
Oxy-degradable			
Bio based Polyethylene			
Starch based - potato			
Starch based - corn			
Starch based - rice			
Starch based - other			
Unwrapped (no outer)			
Paper wrapped			
Paper envelope			

3. Provide your contact details:

Company Name	
Mailmark Participant ID	
Contact Name	
Telephone Number	
Email Address	

4. Send your sample items to the following address:

Mailmark Quality Assurance
Royal Mail
Dearne House
Cortonwood Drive
Brampton
BARNLEY
S73 0UF