

## Large Letter Mailmark barcode quality and item design check

- 1. Provide 10 fully addressed, filled and sealed mail pieces with the appropriate indicia and Mailmark barcode. Items will be retained for future reference
- 2. Provide details of the type of wrapping used for the samples sent.

Please provide a separate sheet per type to enable us to differentiate during testing.

Wrap type	Quantity	Haze	Micron thickness
Polyethylene [PE]			
Polypropylene [PP]			
Oxy-degradable			
Bio based Polyethylene			
Starch based - potato			
Starch based - corn			
Starch based - rice			
Starch based - other			
Unwrapped (no outer)			
Paper wrapped			
Paper envelope			

3. Provide your contact details:

Company Name	
Mailmark Participant ID	
Contact Name	
Telephone Number	
Email Address	

4. Send your sample items to the following address:

Mailmark Quality Assurance Drake House Breakwater Road Plymstock PLYMOUTH PL9 7HW